

Anti-CD83 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0303-FITC

Aliases

CD83; CD83 Molecule; HB15; BL11; CD83 Antigen (Activated B Lymphocytes, Immunoglobulin Superfamily); B-Cell Activation Protein; Cell Surface Protein HB15; CD83 Antigen; HCD83; Cell-Surface Glycoprotein; CD83 Molecule

Background

Gene Name: CD83

NCBI Gene Entry: [9308](#)

UniProt Entry: [Q01151](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14775

Species Reactivity: Human, african green, baboon, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD83 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM