

Anti-CD99 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0307-FITC

Aliases

CD99; CD99 Molecule (Xg Blood Group); CD99 Antigen; MIC2X; MIC2Y; MIC2; Antigen Identified By Monoclonal Antibodies 12E7, F21 And O13; T-Cell Surface Glycoprotein E2; E2 Antigen; Antigen Identified By Monoclonal 12E7, Y Homolog; MIC2 (Monoclonal Antibody 12E7); Cell Surface Antigen HBA-71; Cell Surface Antigen 12E7; Cell Surface Antigen O13; Surface Antigen MIC2; CD99 Molecule; Protein MIC2; HBA71; MSK5X; 12E7

Background

Gene Name: CD99

NCBI Gene Entry: [4267](#)

UniProt Entry: [P14209](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14800

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD99 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM