

Anti-CD103 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0308-FITC

Aliases

ITGAE; Integrin Subunit Alpha E; HUMINAE; CD103; Integrin, Alpha E (Antigen CD103, Human Mucosal Lymphocyte Antigen 1; Alpha Polypeptide); Human Mucosal Lymphocyte Antigen 1, Alpha Polypeptide; Mucosal Lymphocyte 1 Antigen; Integrin Alpha-IEL; Integrin Alpha-E; HML-1 Antigen; Antigen CD103, Human Mucosal Lymphocyte Antigen 1; Alpha Polypeptid; Antigen CD103; CD103 Antigen

Background

Gene Name: ITGAE

NCBI Gene Entry: [3682](#)

UniProt Entry: [P38570](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14805

Species Reactivity: Human, macaque

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD103 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM