

Anti-CD119 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0315-FITC

Aliases

IFNGR1; Interferon Gamma Receptor 1; CD119; IFN-Gamma Receptor 1; IFN-Gamma-R-Alpha; CD119 Antigen; IFN-Gamma-R1; CDw119; IFNGR; Interferon-Gamma Receptor Alpha Chain; Interferon Gamma Receptor Alpha-Chain; Immune Interferon Receptor 1; Antiviral Protein, Type 2; AVP, Type 2; IMD27A; IMD27B

Background

Gene Name: IFNGR1

NCBI Gene Entry: [2459](#)

UniProt Entry: [P15260](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15290

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, k

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

