

# Anti-CD178 Mouse Monoclonal Antibody, FITC Conjugated



**Catalog #: F0325-FITC**

## Aliases

FASLG; Fas Ligand; FasL; APT1LG1; TNFSF6; CD178; Tumor Necrosis Factor Ligand Superfamily Member 6; Fas Ligand (TNF Superfamily, Member 6); Apoptosis Antigen Ligand; Fas Antigen Ligand; CD95 Ligand; CD95-L; CD95L; APTL; Tumor Necrosis Factor (Ligand) Superfamily, Member 6; Mutant Tumor Necrosis Factor Family Member 6; Apoptosis (APO-1) Antigen Ligand 1; Tumor Necrosis Factor Ligand 1A; CD178 Antigen; ALPS1B; TNLG1A; FASL

## Background

Gene Name: FASLG

NCBI Gene Entry: [356](#)

UniProt Entry: [P48023](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14890

Species Reactivity: Human, baboon

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD178 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)