

Anti-CD200 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0326-FITC

Aliases

CD200; CD200 Molecule; CD200 Antigen; OX-2; MOX1; MOX2; MRC; Antigen Identified By Monoclonal Antibody MRC OX-2; OX-2 Membrane Glycoprotein

Background

Gene Name: CD200

NCBI Gene Entry: [4345](#)

UniProt Entry: [P41217](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14905

Species Reactivity: Human, african green, baboon

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD200 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM