

Anti-CD235 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0328-FITC

Aliases

GYPA; Glycophorin A (MNS Blood Group); PAS-2; GPA; Sialoglycoprotein Alpha; MN Sialoglycoprotein; CD235a; MN; Glycophorin A (MN Blood Group); Glycophorin-A; MNS; Erythroid-Lineage-Specific Membrane Sialoglycoprotein; Recombinant Glycophorin A-B Miltenberger-DR; Glycophorin A (Includes MN Blood Group); Mi.V Glycoprotein (24 AA); Glycophorin A Transcript; Glycophorin Sta Type C; Glycophorin A, GPA; Glycophorin Erik; Glycophorin MiI; Glycophorin MiV; Glycophorin SAT; CD235a Antigen; GlycophorinA; HGpSta(C); HGpMiXI; GPErik; HGpMiV; GPSAT

Background

Gene Name: GYPA

NCBI Gene Entry: [2993](#)

UniProt Entry: [P02724](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14915

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD235 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM