

Anti-CD278 Armenian Hamster Monoclonal Antibody, FITC Conjugated



Catalog #: F0330-FITC

Aliases

ICOS; Inducible T Cell Costimulator; AILIM; Activation-Inducible Lymphocyte Immunomediatory Molecule; Inducible T-Cell Costimulator; CD278; Inducible T-Cell Co-Stimulator; Inducible Costimulator; CD278 Antigen; CVID1

Background

Gene Name: ICOS

NCBI Gene Entry: [29851](#)

UniProt Entry: [Q9Y6W8](#)

Application Information

Clonality: Armenian Hamster monoclonal antibody

Clone ID: 24GB14935

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Armenian Hamster IgG

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD278 Armenian Hamster Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM