

Anti-CD303 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0333-FITC

Aliases

CLEC4C; C-Type Lectin Domain Family 4 Member C; BDCA2; HECL; DLEC; CLECSF11; CLECSF7; CD303; C-Type (Calcium Dependent, Carbohydrate-Recognition Domain) Lectin, Superfamily Member 7; C-Type Lectin Superfamily Member 7; Dendritic Lectin; BDCA-2; C-Type (Calcium Dependent, Carbohydrate-Recognition Domain) Lectin, Superfamily Member 11; C-Type Lectin Domain Family 4, Member C; Blood Dendritic Cell Antigen 2 Protein; Blood Dendritic Cell Antigen 2; Dendritic Cell Lectin B; CD303 Antigen; PRO34150

Background

Gene Name: CLEC4C

NCBI Gene Entry: [170482](#)

UniProt Entry: [Q8WTT0](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14945

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD303 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM