

Anti-CD314 Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0334-FITC

Aliases

TNFRSF4; TNF Receptor Superfamily Member 4; TXGP1L; ACT35; CD134; OX40; TAX Transcriptionally-Activated Glycoprotein 1 Receptor; Tumor Necrosis Factor Receptor Superfamily Member 4; OX40L Receptor; ACT35 Antigen; CD134 Antigen; Tax-Transcriptionally Activated Glycoprotein 1 Receptor; Tumor Necrosis Factor Receptor Superfamily, Member 4; Lymphoid Activation Antigene ACT35; OX40 Cell Surface Antigen; OX40 Homologue; ATC35 Antigen; OX40 Antigen; IMD16

Background

Gene Name: TNFRSF4

NCBI Gene Entry: [7293](#)

UniProt Entry: [P43489](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14955

Species Reactivity: Human, african green, baboon, cynomolgus, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD314 Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM