

Anti-HLA-DP Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0338-FITC

Aliases

HLA-DPA1; Major Histocompatibility Complex, Class II, DP Alpha 1; HLA-DP1A; HLA Class II Histocompatibility Antigen, DP Alpha 1 Chain; MHC Class II DP3-Alpha; HLA-SB Alpha Chain; DP(W3); DP(W4); HLASB; MHC Class II Antigen DP Alpha Chain; MHC Class II HLA-DPA1 Antigen; MHC Class II DP Alpha Chain; MHC Class II Antigen; MHC Class II Protein; MHC Class II DPA1; HLA-DPB1; HLA DPA1; HLA-DPA; HLADP; DPA1; PLT1; HLA-DPB1; Major Histocompatibility Complex, Class II, DP Beta 1; HLA-DP1B; HLA Class II Histocompatibility Antigen, DP(W4) Beta Chain; HLA Class II Histocompatibility Antigen, DP Beta 1 Chain; MHC Class II Antigen DPB1; HLA-DP Histocompatibility Type, Beta-1 Subunit; MHC Class II HLA-DP-Beta-1; MHC HLA DPB1; HLA-DPB; HLA-DP; DPB1

Background

Gene Name: HLA-DPA1/HLA-DPB1

NCBI Gene Entry: [3113/3115](#)

UniProt Entry: [P20036/P04440](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15325

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG3

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM


Anti-HLA-DP Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2

5 µL/Tests

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM