

# Anti-HLA-DQ Mouse Monoclonal Antibody, FITC Conjugated



## Catalog #: F0339-FITC

### Aliases

HLA-DQA1; Major Histocompatibility Complex, Class II, DQ Alpha 1; CELIAC1; HLA Class II Histocompatibility Antigen, DQ Alpha 1 Chain; MHC Class II DQA1; DC-1 Alpha Chain; DC-Alpha; HLA-DQA; HLA-DCA; HLA Class II Histocompatibility Antigen DQ Alpha Chain; MHC Class II DQ Alpha Chain; MHC Class II HLA-DQ-Alpha-1; MHC Class II Antigen DQA1; MHC Class II Protein; MHC HLA-DQ Alpha; HLA-DQA1\*; HLA-DQB1; DQ-A1; DQA1; HLA-DQB1; Major Histocompatibility Complex, Class II, DQ Beta 1; CELIAC1; HLA-DQB; IDDM1; HLA Class II Histocompatibility Antigen, DQ Beta 1 Chain; MHC Class II Antigen DQB1; MHC Class II HLA-DQ Beta Glycoprotein; MHC Class II Antigen HLA-DQ-Beta-1; MHC Class II DQ Beta Chain

### Background

Gene Name: HLA-DQA1/HLA-DQB1

NCBI Gene Entry: [3117/3119](#)

UniProt Entry: [P01909/P01920](#)

### Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15330

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

### Immunogen

The immunogen information is proprietary.

### Isotype

Mouse IgG2a

### Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

### Storage

Store at 4 °C and protected from prolonged exposure to light.

### Recommended Dilutions

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX:+1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-HLA-DQ Mouse Monoclonal Antibody, FITC Conjugated

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data



---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)