

Anti-IFN- γ Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0343-FITC

Aliases

IFNG; Interferon Gamma; Immune Interferon; IFN-Gamma; IMD69; IFG; IFI

Background

Gene Name: IFNG

NCBI Gene Entry: [3458](#)

UniProt Entry: [P01579](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14970

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

