

Anti-IL1-beta Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0344-FITC

Aliases

IL1B; Interleukin 1 Beta; IL1F2; IL1-BETA; Interleukin-1 Beta; IL-1 Beta; Catabolin; IL-1B; Pro-Interleukin-1-Beta; Preinterleukin 1 Beta; Interleukin 1beta; IL1beta; IL-1

Background

Gene Name: IL1B

NCBI Gene Entry: [3553](#)

UniProt Entry: [P01584](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15345

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-IL1-beta Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM