

Anti-IL-12p40 Monoclonal Antibody, FITC Conjugated



Catalog #: F0345-FITC

Aliases

IL12B; Interleukin 12B; IL-12B; CLMF2; NKSF2; CLMF; NKSF; Interleukin 12B (Natural Killer Cell Stimulatory Factor 2, Cytotoxic Lymphocyte Maturation Factor 2, P40); Natural Killer Cell Stimulatory Factor, 40 KD Subunit; Cytotoxic Lymphocyte Maturation Factor 40 KDa Subunit; NK Cell Stimulatory Factor Chain 2; Interleukin-12 Subunit Beta; Interleukin-12 Beta Chain; Interleukin 12, P40; IL12, Subunit P40; IL-12 Subunit P40; CLMF P40; Cytotoxic Lymphocyte Maturation Factor 2, P40; Natural Killer Cell Stimulatory Factor-2; IMD28; IMD29

Background

Gene Name: IL12B

NCBI Gene Entry: [3593](#)

UniProt Entry: [P29460](#)

Application Information

Clonality: Monoclonal antibody

Clone ID: 24GB14980

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Human IgG1, λ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-IL-12p40 Monoclonal Antibody, FITC Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM