

Anti-Lambda Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0350-FITC

Aliases

POLL; DNA Polymerase Lambda; Polymerase (DNA Directed), Lambda; Polymerase (DNA) Lambda; DNA Polymerase Beta-2; DNA Polymerase Kappa; DNA Polymerase Beta-N; EC 4.2.99.-; Pol Lambda; EC 2.7.7.7; Pol Beta2; POLKAPPA; BETAN

Background

Gene Name: POLL

NCBI Gene Entry: [27343](#)

UniProt Entry: [Q9UGP5](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15365

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Lambda Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM