

Anti-myeloperoxidase Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0351-FITC

Aliases

MPO; Myeloperoxidase; EC 1.11.2.2; EC 1.11.1.7; EC 1.11.1

Background

Gene Name: MPO

NCBI Gene Entry: [4353](#)

UniProt Entry: [P05164](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15370

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-myeloperoxidase Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM