

Anti-TCR γ/δ Mouse Monoclonal Antibody, FITC Conjugated



Catalog #: F0353-FITC

Aliases

TRD; T Cell Receptor Delta Locus; TCRD, TRD@; TCRDV1; T-Cell Antigen Receptor, Delta Polypeptide; T-Cell Receptor, Delta (V,D,J,C); TCRD; TRD@; TRG; T Cell Receptor Gamma Locus; T-Cell Antigen Receptor, Gamma Polypeptide; T-Cell Receptor, Gamma Cluster; T-Cell Rearranging Gene, Gamma; TCRG, TRG@; TCRG; TRG@

Background

Gene Name: TRD/TRG

NCBI Gene Entry: [6964/6965](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15375

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-TCR γ/δ Mouse Monoclonal Antibody, FITC Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM