

Anti-CD11c Mouse Monoclonal Antibody, AF488 Conjugated



Catalog #: F0358-AF488

Aliases

ITGAX; Integrin Subunit Alpha X; CD11c; Integrin, Alpha X (Complement Component 3 Receptor 4 Subunit); Integrin, Alpha X (Antigen CD11C (P150), Alpha Polypeptide); Leukocyte Adhesion Glycoprotein P150,95 Alpha Chain; Complement Component 3 Receptor 4 Subunit; Leukocyte Adhesion Receptor P150,95; CD11 Antigen-Like Family Member C; Integrin Alpha-X; Leukocyte Surface Antigen P150,95, Alpha Subunit; Myeloid Membrane Antigen, Alpha Subunit; P150 95 Integrin Alpha Chain; Leu M5, Alpha Subunit; Integrin Alpha X; CD11c Antigen; Leu M5; SLEB6; CD11C

Background

Gene Name: ITGAX

NCBI Gene Entry: [3687](#)

UniProt Entry: [P20702](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14475

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: AF488

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD11c Mouse Monoclonal Antibody, AF488 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM