

Anti-Ki-67 Mouse Monoclonal Antibody, AF488 Conjugated



Catalog #: F0369-AF488

Aliases

MKI67; Marker Of Proliferation Ki-67; Antigen Identified By Monoclonal Antibody Ki-67; PPP1R105; MIB-1; Protein Phosphatase 1, Regulatory Subunit 105; Molecular Immunology Borstel Antibody 1; Proliferation Marker Protein Ki-67; Antigen Ki67; Ki-67; Proliferation-Related Ki-67 Antigen; Antigen KI-67; MIB-; KIA

Background

Gene Name: MKI67

NCBI Gene Entry: [4288](#)

UniProt Entry: [P46013](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15400

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: AF488

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Ki-67 Mouse Monoclonal Antibody, AF488 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM