

# Anti-CD11c Mouse Monoclonal Antibody, AF700 Conjugated



**Catalog #: F0370-AF700**

## Aliases

ITGAX; Integrin Subunit Alpha X; CD11c; Integrin, Alpha X (Complement Component 3 Receptor 4 Subunit); Integrin, Alpha X (Antigen CD11C (P150), Alpha Polypeptide); Leukocyte Adhesion Glycoprotein P150,95 Alpha Chain; omplement Component 3 Receptor 4 Subunit; Leukocyte Adhesion Receptor P150,95; CD11 Antigen-Like Family Member C; Integrin Alpha-X; Leukocyte Surface Antigen P150,95, Alpha Subunit; Myeloid Membrane Antigen, Alpha Subunit; P150 95 Integrin Alpha Chain; Leu M5, Alpha Subunit; Integrin Alpha X; CD11c Antigen; Leu M5; SLEB6; CD11C

## Background

Gene Name: ITGAX

NCBI Gene Entry: [3687](#)

UniProt Entry: [P20702](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14475

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: AF700

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD11c Mouse Monoclonal Antibody, AF700 Conjugated

PAGE 2

**Note:** This product is for research use only.

### Validation Data



---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)