

# Anti-CD45 Mouse Monoclonal Antibody, AF647 Conjugated



**Catalog #: F0387-AF647**

## Aliases

PTPRC; Protein Tyrosine Phosphatase Receptor Type C; T200; GP180; CD45R; B220; CD45; LCA; LY5; Receptor-Type Tyrosine-Protein Phosphatase C; Glycoprotein 180; CD45 Antigen; LCA; Protein Tyrosine Phosphatase, Receptor Type, C Polypeptide; T200 Leukocyte Common Antigen; Leukocyte-Common Antigen; Leukocyte Common Antigen; T200 Glycoprotein; EC 3.1.3.48; IMD105

## Background

Gene Name: PTPRC  
NCBI Gene Entry: [5788](#)  
UniProt Entry: [P08575](#)

## Application Information

Clonality: Mouse monoclonal antibody  
Clone ID: 24GB14650  
Species Reactivity: Human  
Applications Tested: Flow Cytometry (FCM)  
Conjugation: AF647

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG2a

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD45 Mouse Monoclonal Antibody, AF647 Conjugated

## Validation Data

