

Anti-CD158e1 Mouse Monoclonal Antibody, AF647 Conjugated



Catalog #: F0392-AF647

Aliases

KIR3DL1; Killer Cell Immunoglobulin Like Receptor, Three Ig Domains And Long Cytoplasmic Tail 1; NKB1; CD158E1; Nkat3; NKB1B; Killer Cell Immunoglobulin-Like Receptor, Three Domains, Long Cytoplasmic Tail, 1; P70 Natural Killer Cell Receptor Clones CL-2/CL-11; Killer Cell Immunoglobulin-Like Receptor 3DL1; HLA-BW4-Specific Inhibitory NK Cell Receptor; Natural Killer-Associated Transcript 3; CD158 Antigen-Like Family Member E; P70 NK Receptor CL-2/CL-11; CD158e1/2; CD158e2; NKAT-3; Cl-11; AMB11; Cl-2; KIR; P70 Killer Cell Inhibitory Receptor; KIR Antigen 3DL1; CD158e Antigen; KIR3DL1/S1; KIR2DL5B; CD158E; NKAT3

Background

Gene Name: KIR3DL1

NCBI Gene Entry: [3811](#)

UniProt Entry: [P43629](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14880

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: AF647

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD158e1 Mouse Monoclonal Antibody, AF647 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM