

Anti-Ki-67 Mouse Monoclonal Antibody, AF647 Conjugated



Catalog #: F0394-AF647

Aliases

MKI67; Marker Of Proliferation Ki-67; Antigen Identified By Monoclonal Antibody Ki-67; PPP1R105; MIB-1; Protein Phosphatase 1, Regulatory Subunit 105; Molecular Immunology Borstel Antibody 1; Proliferation Marker Protein Ki-67; Antigen Ki67; Ki-67; Proliferation-Related Ki-67 Antigen; Antigen KI-67; MIB-; KIA

Background

Gene Name: MKI67
NCBI Gene Entry: 4288
UniProt Entry: P46013

Application Information

Clonality: Mouse monoclonal antibody
Clone ID: 24GB15400
Species Reactivity: Human
Applications Tested: Intracellular Staining for Flow Cytometry ICFC
Conjugation: AF647

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT


SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Ki-67 Mouse Monoclonal Antibody, AF647 Conjugated



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM