

# Anti-CD4 Mouse Monoclonal Antibody, PE-Cy7 Conjugated



**Catalog #: F0398-PC7**

## Aliases

CD4 □ CD4 Molecule □ T-Cell Surface Glycoprotein CD4 □ Leu-3 □ T4 □ T-Cell Surface Antigen T4/Leu-3 □ CD4 Antigen (P55) □ CD4 Receptor □ CD4 Antigen □ CD4mut □ IMD79 □ OKT4D

## Background

Gene Name: CD4

NCBI Gene Entry: [920](#)

UniProt Entry: [P01730](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14420

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG2b

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT


SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD4 Mouse Monoclonal Antibody, PE-Cy7 Conjugated



**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)