

Anti-CD7 Mouse Monoclonal Antibody, PE-Cy7 Conjugated



Catalog #: F0401-PC7

Aliases

CD7; CD7 Molecule; GP40; TP41; T-Cell Leukemia Antigen; T-Cell Antigen CD7; LEU-9; Tp40; T-Cell Surface Antigen Leu-9 ; CD7 Antigen (P41); P41 Protein; CD7 Antigen

Background

Gene Name: CD7

NCBI Gene Entry: [924](#)

UniProt Entry: [P09564](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14435

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD7 Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM