

Anti-CD11b Mouse Monoclonal Antibody, PE-Cy7 Conjugated



Catalog #: F0409-PC7

Aliases

ITGAM; Integrin Subunit Alpha M; MAC-1; CD11b; CR3A; Integrin, Alpha M (Complement Component 3 Receptor 3 Subunit); Cell Surface Glycoprotein MAC-1 Subunit Alpha; Complement Component 3 Receptor 3 Subunit; Macrophage-1 Antigen Alpha Subunit; CD11 Antigen-Like Family Member B; Leukocyte Adhesion Receptor MO1; Integrin Alpha-M; CR-3 Alpha Chain; Integrin, Alpha M (Complement Component Receptor 3, Alpha; Also Known As CD11b (P170), Macrophage Antigen Alpha Polypeptide); Neutrophil Adherence Receptor Alpha-M Subunit; Macrophage Antigen Alpha Polypeptide; Neutrophil Adherence Receptor; Antigen CD11b (P170); CD11b Antigen; MAC1A; SLEB6; CD11B; MO1A

Background

Gene Name: ITGAM

NCBI Gene Entry: [3684](#)

UniProt Entry: [P11215](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15425

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD11b Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM