

Anti-CD15 Mouse Monoclonal Antibody, PE-Cy7 Conjugated



Catalog #: F0411-PC7

Aliases

FUT4; Fucosyltransferase 4; Galactoside 3-L-Fucosyltransferase; FUC-TIV; FCT3A; ELFT; 4-Galactosyl-N-Acetylglucosaminide 3-Alpha-L-Fucosyltransferase; Alpha (1,3) Fucosyltransferase, Myeloid-Specific; Alpha-(1,3)-Fucosyltransferase 4; ELAM-1 Ligand Fucosyltransferase; ELAM Ligand Fucosyltransferase; Fucosyltransferase IV; FucT-IV; CD15; Stage-Specific Embryonic Antigen 1; EC 2.4.1.152; EC 2.4.1.65; EC 2.4.1; Lewis X; Fuc-TIV; SSEA-1; FUTIV; LeX

Background

Gene Name: FUT4

NCBI Gene Entry: [2526](#)

UniProt Entry: [P22083](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14495

Species Reactivity: Human, chimpanzee

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgM

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD15 Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM