

Anti-CD23 Mouse Monoclonal Antibody, PE-Cy7 Conjugated



Catalog #: F0416-PC7

Aliases

FCER2; Fc Epsilon Receptor II; CLEC4J; FcepsilonRII; FCErII; CD23A; CD23; FCE2; Fc Fragment Of IgE, Low Affinity II, Receptor For (CD23); Low Affinity Immunoglobulin Epsilon Fc Receptor; C-Type Lectin Domain Family 4 Member J; Immunoglobulin E-Binding Factor; Fc Fragment Of IgE Receptor II; Lymphocyte IgE Receptor; Fc-Epsilon-RII; CD23 Antigen; BLAST-2; IGEBF; Fc Fragment Of IgE, Low Affinity II, Receptor For (CD23A); C-Type Lectin Domain Family 4, Member J; Immunoglobulin Epsilon-Chain; Fc Epsilon RII

Background

Gene Name: FCER2

NCBI Gene Entry: [2208](#)

UniProt Entry: [P06734](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14540

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD23 Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM