

# Anti-CD25 Mouse Monoclonal Antibody, PE-Cy7 Conjugated



**Catalog #: F0418-PC7**

## Aliases

IL2RA; Interleukin 2 Receptor Subunit Alpha; CD25; Interleukin-2 Receptor Subunit Alpha; Interleukin 2 Receptor, Alpha; IL-2 Receptor Subunit Alpha; IL-2R Subunit Alpha; TAC Antigen; IDDM10; IL2R; P55; Insulin-Dependent Diabetes Mellitus 10; CD25 Antigen; IL-2-RA; IL2-RA; IMD41; TCGFR

## Background

Gene Name: IL2RA

NCBI Gene Entry: [3559](#)

UniProt Entry: [P01589](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15385

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-CD25 Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)