

Anti-CD29 Mouse Monoclonal Antibody, PE-Cy7 Conjugated



Catalog #: F0419-PC7

Aliases

ITGB1; Integrin Subunit Beta 1; GPIIA; MSK12; CD29; FNRB; MDF2; Integrin, Beta 1 (Fibronectin Receptor, Beta Polypeptide, Antigen CD29 Includes MDF2, MSK12); Glycoprotein IIa; Integrin Beta-1; Very Late Activation Protein, Beta Polypeptide; Fibronectin Receptor Subunit Beta; Integrin VLA-4 Beta Subunit; VLA-4 Subunit Beta; Integrin Beta 1; CD29 Antigen; VLA-BETA; VLAB

Background

Gene Name: ITGB1

NCBI Gene Entry: [3688](#)

UniProt Entry: [P05556](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14580

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD29 Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM