

Anti-CD79a Mouse Monoclonal Antibody, PE-Cy7 Conjugated



Catalog #: F0431-PC7

Aliases

CD79A; CD79a Molecule; Ig-Alpha; MB1; B-Cell Antigen Receptor Complex-Associated Protein Alpha Chain; IGAAlpha; MB-1; IGA; Membrane-Bound Immunoglobulin-Associated Protein; CD79A Antigen (Immunoglobulin-Associated Alpha); CD79a Molecule, Immunoglobulin-Associated Alpha; Surface IgM-Associated Protein; MB-1 Membrane Glycoprotein; CD79a Antigen

Background

Gene Name: CD79A

NCBI Gene Entry: [973](#)

UniProt Entry: [P11912](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14765

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD79a Mouse Monoclonal Antibody, PE- Cy7 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM