

Anti-CD79b Mouse Monoclonal Antibody, PE-Cy7 Conjugated



Catalog #: F0432-PC7

Aliases

CD79B; CD79b Molecule; Ig-Beta; B29; B-Cell Antigen Receptor Complex-Associated Protein Beta Chain; Igbeta; IGB; CD79B Antigen (Immunoglobulin-Associated Beta); CD79b Molecule, Immunoglobulin-Associated Beta; Immunoglobulin-Associated B29 Protein; B-Cell-Specific Glycoprotein B29; CD79b Antigen; AGM6

Background

Gene Name: CD79B

NCBI Gene Entry: [974](#)

UniProt Entry: [P40259](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15090

Species Reactivity: Human, african green, baboon, cynomolgus, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD79b Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM