

# Anti-TCR $\gamma/\delta$ Mouse Monoclonal Antibody, PE-Cy7 Conjugated



**Catalog #: F0441-PC7**

## Aliases

TRD; T Cell Receptor Delta Locus; TCRD, TRD@; TCRDV1; T-Cell Antigen Receptor, Delta Polypeptide; T-Cell Receptor, Delta (V,D,J,C); TCRD; TRD@; TRG; T Cell Receptor Gamma Locus; T-Cell Antigen Receptor, Gamma Polypeptide; T-Cell Receptor, Gamma Cluster; T-Cell Rearranging Gene, Gamma; TCRG, TRG@; TCRG; TRG@

## Background

Gene Name: TRD/TRG

NCBI Gene Entry: [6964/6965](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15375

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE-Cy7

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-TCR $\gamma/\delta$ Mouse Monoclonal Antibody, PE-Cy7 Conjugated

PAGE 2



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)