

Anti-CD13 Mouse Monoclonal Antibody, APC-Cy7 Conjugated



Catalog #: F0451-AC7

Aliases

ANPEP; Alanyl Aminopeptidase, Membrane; Gp150; AP-N; HAPN; Microsomal Aminopeptidase; Aminopeptidase N; Aminopeptidase M; LAP1; P150; CD13; PEPN; Myeloid Plasma Membrane Glycoprotein CD13; Alanyl (Membrane) Aminopeptidase; Membrane Alanyl Aminopeptidase; EC 3.4.11.2; AP-M; APN; Alanyl Aminopeptidase; CD13 Antigen; EC 3.4.11

Background

Gene Name: ANPEP

NCBI Gene Entry: [290](#)

UniProt Entry: [P15144](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14480

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD13 Mouse Monoclonal Antibody, APC-Cy7 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM