

Anti-CD36 Mouse Monoclonal Antibody, APC-Cy7 Conjugated



Catalog #: F0457-AC7

Aliases

CD36; CD36 Molecule (CD36 Blood Group); GPIIIB; GPIV; GP3B; FAT; GP4; Platelet Glycoprotein 4; Fatty Acid Translocase; SCARB3; CD36 Antigen (Collagen Type I Receptor, Thrombospondin Receptor); CD36 Molecule (Thrombospondin Receptor); Leukocyte Differentiation Antigen CD36; Platelet Glycoprotein IV; Glycoprotein IIb; PAS IV; Scavenger Receptor Class B, Member 3; Scavenger Receptor Class B Member 3; Platelet Collagen Receptor; Thrombospondin Receptor; Cluster Determinant 36; PAS-4 Protein; CD36 Antigen; BDPLT10; CHDS7; PASIV; PAS-4

Background

Gene Name: CD36

NCBI Gene Entry: [948](#)

UniProt Entry: [P16671](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14610

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD36 Mouse Monoclonal Antibody, APC-Cy7 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM