

Anti-CD41a Mouse Monoclonal Antibody, APC-Cy7 Conjugated



Catalog #: F0459-AC7

Aliases

ITGA2B; Integrin Subunit Alpha 2b; GPIIb; Integrin Alpha-IIb; PPP1R93; CD41B; CD41; GP2B; Integrin, Alpha 2b (Platelet Glycoprotein IIb Of IIb/IIIa Complex, Antigen CD41); Platelet Glycoprotein IIb Of IIb/IIIa Complex; Protein Phosphatase 1, Regulatory Subunit 93; Platelet Membrane Glycoprotein IIb; GPalpha IIb; Integrin, Alpha 2b (Platelet Glycoprotein IIb Of IIb/IIIa Complex, Antigen CD41B); Platelet Fibrinogen Receptor, Alpha Subunit; Platelet-Specific Antigen BAK; AlphaIIb Protein; CD41 Antigen; BDPLT16; BDPLT2; ITGAB; HPA3; GT1; GTA; GT

Background

Gene Name: ITGA2B

NCBI Gene Entry: [3674](#)

UniProt Entry: [P08514](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14625

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD41a Mouse Monoclonal Antibody, APC-Cy7 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM