

Anti-CD45RO Mouse Monoclonal Antibody, APC-Cy7 Conjugated



Catalog #: F0462-AC7

Aliases

PTPRC; Protein Tyrosine Phosphatase Receptor Type C; T200; GP180; CD45R; B220; CD45; LCA; LY5; Receptor-Type Tyrosine-Protein Phosphatase C; Glycoprotein 180; CD45 Antigen; LCA; Protein Tyrosine Phosphatase, Receptor Type, C Polypeptide; T200 Leukocyte Common Antigen; Leukocyte-Common Antigen; Leukocyte Common Antigen; T200 Glycoprotein; EC 3.1.3.48; IMD105

Background

Gene Name: PTPRC

NCBI Gene Entry: [5788](#)

UniProt Entry: [P08575](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14660

Species Reactivity: Human, chimpanzee, cynomolgus, marmoset

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD45RO Mouse Monoclonal Antibody, APC-Cy7 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM