

Anti-CD49d Mouse Monoclonal Antibody, APC-Cy7 Conjugated



Catalog #: F0463-AC7

Aliases

ITGA4; Integrin Subunit Alpha 4; CD49d; Alpha 4 Subunit Of VLA-4 Receptor; CD49 Antigen-Like Family Member D; VLA-4 Subunit Alpha; Integrin Alpha-IV; Integrin Alpha-4; Integrin, Alpha 4 (Antigen CD49D, Alpha 4 Subunit Of VLA-4 Receptor); Very Late Activation Protein 4 Receptor, Alpha 4 Subunit; Antigen CD49D, Alpha-4 Subunit Of VLA-4 Receptor; 269C Wild Type; Antigen CD49D; CD49d Antigen; CD49D; IA4

Background

Gene Name: ITGA4

NCBI Gene Entry: [3676](#)

UniProt Entry: [P13612](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15215

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD49d Mouse Monoclonal Antibody, APC-Cy7 Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM