

Anti-HLA-DR Mouse Monoclonal Antibody, APC-Cy7 Conjugated



Catalog #: F0466-AC7

Aliases

HLA-DRA; Major Histocompatibility Complex, Class II, DR Alpha; HLA-DRA1; HLA Class II Histocompatibility Antigen, DR Alpha Chain; MHC Class II Antigen DRA; Histocompatibility Antigen HLA-DR Alpha; HLA-DRB1; Major Histocompatibility Complex, Class II, DR Beta 1; HLA-DR1B; Major Histocompatibility Complex, Class II, DR Beta 1 Precursor; HLA Class II Histocompatibility Antigen, DR-1 Beta Chain; HLA Class II Histocompatibility Antigen, DRB1 Beta Chain; MHC Class II HLA-DR Beta 1 Chain; Human Leucocyte Antigen DRB1; Human Leukocyte Antigen DRB1; Lymphocyte Antigen DRB1; HLA-DRB; DRB1; SS1

Background

Gene Name: HLA-DRA/HLA-DRB1

NCBI Gene Entry: [3122/3123](#)

UniProt Entry: [P01903/P01911](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14965

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-HLA-DR Mouse Monoclonal Antibody, APC-Cy7 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM