

# Anti-CD3 Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0470-PE**

## Aliases

CD3G; CD3 Gamma Subunit Of T-Cell Receptor Complex; T-Cell Surface Glycoprotein CD3 Gamma Chain; CD3-GAMMA; CD3GAMMA; CD3g Antigen, Gamma Polypeptide (TiT3 Complex); CD3g Molecule, Gamma (CD3-TCR Complex); T-Cell Receptor T3 Gamma Chain; T3G; T-Cell Antigen Receptor Complex, Gamma Subunit Of T3; CD3g Molecule, Epsilon (CD3-TCR Complex); CD3g Molecule; CD3g Antigen; IMD17

## Background

Gene Name: CD3G

NCBI Gene Entry: [917](#)

UniProt Entry: [P09693](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14405

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG2a

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

**Note:** This product is for research use only.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD3 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)