

# Anti-CD18 Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0490-PE**

## Aliases

ITGB2; Integrin Subunit Beta 2; LFA-1; MAC-1; CD18; MFI7; Integrin, Beta 2 (Complement Component 3 Receptor 3 And 4 Subunit); Complement Component 3 Receptor 3 And 4 Subunit; Integrin Beta-2; Integrin, Beta 2 (Antigen CD18 (P95), Lymphocyte Function-Associated Antigen 1; Macrophage Antigen 1 (Mac-1) Beta Subunit); Cell Surface Adhesion Glycoprotein (LFA-1/CR3/P150,959 Beta Subunit Precursor); Cell Surface Adhesion Glycoproteins LFA-1/CR3/P150,95 Subunit Beta; Leukocyte-Associated Antigens CD18/11A, CD18/11B, CD18/11C; Leukocyte Cell Adhesion Molecule CD18; Complement Receptor C3 Beta-Subunit; Complement Receptor C3 Subunit Beta; Integrin Beta Chain, Beta 2; D18 Antigen; LCAMB; MF17; LAD

## Background

Gene Name: ITGB2

NCBI Gene Entry: [3689](#)

UniProt Entry: [P05107](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14505

Species Reactivity: Human, pig

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD18 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)