

Anti-CD40 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0517-PE

Aliases

CD40; CD40 Molecule; Bp50; TNFRSF5; P50; Tumor Necrosis Factor Receptor Superfamily Member 5; CD40 Molecule, TNF Receptor Superfamily Member 5; CD40L Receptor; Tumor Necrosis Factor Receptor Superfamily, Member 5; B Cell Surface Antigen CD40; B-Cell Surface Antigen CD40; B Cell-Associated Molecule; CD40 Antigen; CDW40; CDw40

Background

Gene Name: CD40

NCBI Gene Entry: [958](#)

UniProt Entry: [P25942](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14620

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD40 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM