

Anti-CD42a Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0519-PE

Aliases

GP9; Glycoprotein IX Platelet; GPIX; Platelet Glycoprotein IX; CD42a; Glycoprotein 9; Glycoprotein IX (Platelet); CD42a Antigen; GP-IX

Background

Gene Name: GP9

NCBI Gene Entry: [2815](#)

UniProt Entry: [P14770](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14630

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD42a Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM