

# Anti-CD45 Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0525-PE**

## Aliases

PTPRC; Protein Tyrosine Phosphatase Receptor Type C; T200; GP180; CD45R; B220; CD45; LCA; LY5; Receptor-Type Tyrosine-Protein Phosphatase C; Glycoprotein 180; CD45 Antigen; LCA; Protein Tyrosine Phosphatase, Receptor Type, C Polypeptide; T200 Leukocyte Common Antigen; Leukocyte-Common Antigen; Leukocyte Common Antigen; T200 Glycoprotein; EC 3.1.3.48; IMD105

## Background

Gene Name: PTPRC

NCBI Gene Entry: [5788](#)

UniProt Entry: [P08575](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15190

Species Reactivity: Human, chimpanzee

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG1,  $\kappa$

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5  $\mu$ L/Tests

**Note:** This product is for research use only.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD45 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)