

Anti-CD48 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0531-PE

Aliases

CD48; CD48 Molecule; SLAMF2; BLAST; MCD48; HCD48; BCM1; Signaling Lymphocytic Activation Molecule 2; CD48 Antigen (B-Cell Membrane Protein); B-Lymphocyte Activation Marker BLAST-1; Leukocyte Antigen MEM-102; BCM1 Surface Antigen; SLAM Family Member 2; CD48 Antigen; BLAST1; TCT.1; CD48 Molecule; MEM-102

Background

Gene Name: CD48

NCBI Gene Entry: [962](#)

UniProt Entry: [P09326](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15205

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD48 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM