

# Anti-CD49b Mouse Monoclonal Antibody, PE Conjugated



**Catalog #: F0532-PE**

## Aliases

ITGA2; Integrin Subunit Alpha 2; GPIa; VLAA2; HPA-5; CD49B; Integrin, Alpha 2 (CD49B, Alpha 2 Subunit Of VLA-2 Receptor); Very Late Activation Receptor Subunit Alpha-2; Alpha 2 Subunit Of VLA-2 Receptor; CD49 Antigen-Like Family Member B; Platelet Membrane Glycoprotein Ia; Human Platelet Alloantigen 5; Collagen Receptor; Alpha 2 Integrin; Integrin Alpha-2; Glycoprotein Ia; Very Late Activation Protein 2 Receptor, Alpha-2 Subunit; Very Late Activation Receptor Alpha-2 Subunit; Human Platelet Alloantigen System 5; Platelet Glycoprotein GPIa; Platelet Antigen Br; VLA-2 Subunit Alpha; CD49b Antigen; VLA-2; BR

## Background

Gene Name: ITGA2

NCBI Gene Entry: [3673](#)

UniProt Entry: [P17301](#)

## Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15465

Species Reactivity: Human

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

## Immunogen

The immunogen information is proprietary.

## Isotype

Mouse IgG2a

## Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

## Storage

Store at 4 °C and protected from prolonged exposure to light.

## Recommended Dilutions

5 µL/Tests

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

# Anti-CD49b Mouse Monoclonal Antibody, PE Conjugated

PAGE 2

**Note:** This product is for research use only.

## Validation Data



### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)