

Anti-CD68 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0555-PE

Aliases

CD68; CD68 Molecule; CD68 Antigen; Macrosialin; SCARD1; GP110; LAMP4; Scavenger Receptor Class D, Member 1; Macrophage Antigen CD68; DKFZp686M18236; Gp110

Background

Gene Name: CD68

NCBI Gene Entry: [968](#)

UniProt Entry: [P34810](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14745

Species Reactivity: Human

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2b, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD68 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM