

Anti-CD73 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0558-PE

Aliases

NT5E; 5'-Nucleotidase Ecto; CALJA; CD73; ENT; NT5; EN; IMP-Specific 5'-Nucleotidase; Thymidylate 5'-Phosphatase; 5'-Deoxynucleotidase; Ecto-5'-Nucleotidase; 5'-Nucleotidase; EC 3.1.3.5; 5'-NT; NTE; 5'-Nucleotidase, Ecto (CD73); Purine 5-Prime-Nucleotidase; 5' Nucleotidase (CD73); CD73 Antigen; EC 3.1.3.35; EC 3.1.3.89; EC 3.1.3.91; EC 3.1.3.99; E5NT; NT

Background

Gene Name: NT5E

NCBI Gene Entry: [4907](#)

UniProt Entry: [P21589](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14760

Species Reactivity: Human, african green, baboon

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD73 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM