

Anti-CD86 Mouse Monoclonal Antibody, PE Conjugated



Catalog #: F0564-PE

Aliases

CD86; CD86 Molecule; CD28LG2; B7.2; B7-2; CD86 Antigen (CD28 Antigen Ligand 2, B7-2 Antigen); T-Lymphocyte Activation Antigen CD86; CTLA-4 Counter-Receptor B7.2; FUN-1; BU63; B70; B-Lymphocyte Activation Antigen B7-2; B-Lymphocyte Antigen B7-2; Activation B7-2 Antigen; CD86 Antigen; LAB72

Background

Gene Name: CD86

NCBI Gene Entry: [942](#)

UniProt Entry: [P42081](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB14780

Species Reactivity: Human, african green, baboon, rhesus

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD86 Mouse Monoclonal Antibody, PE Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM